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Contents

1	Intended use	27
2	Safety	28
2.1	Safety information	28
2.2	General information	29
2.3	Warnings in this document	30
3	Product description	31
3.1	Module overview	31
3.2	Accessories	32
3.3	Labels and symbols	33
4	Preparation and operation	35
4.1	Connecting the module up to the therapy device	35
4.2	Removing the module from the therapy device	36
4.3	Connecting the pulse oximetry sensor	37
4.4	Connecting the nurse call system and remote alarm box	38
5	Hygienic preparation	40
5.1	General information	40
5.2	Intervals	40
5.3	Hygienic preparation of the module	40
6	Function check	41
6.1	Intervals	41
6.2	Performing a function check	41
7	Faults	42
8	Maintenance	42
9	Storage and disposal	43
9.1	Storage	43
9.2	Disposal	43
10	Appendix	44
10.1	Technical data	44
10.2	Safety distances	45
10.3	Scope of supply	45
10.4	Accessories	45

10.5 Spare parts 46

10.6 Warranty 46

10.7 Declaration of conformity 46

1 Intended use

The WM 110 MS module is used to determine SpO₂ and pulse rate data in combination with compatible therapy devices. A clinical assessment must not be based purely on the values from the module.

The module is also used to connect the therapy device to an intercom system.

WM 100 TD therapy devices do not support connection to an intercom system.

2 Safety

Please read these instructions for use carefully. They form part of the devices described, and must be available at all times.

Use the unit for the designated purpose only (see "1 Intended use", page 27).

For your own safety and that of your patients, and in accordance with the requirements of Directive 93/42/EEC, please observe the following safety instructions.

2.1 Safety information

Warning

Risk of injury due to device or component malfunction!

A damaged device or damaged components may result in injury to the patient, user or bystanders.

⇒ Only operate the device and components if they are externally undamaged.

⇒ Only operate the device and components if the function check has been successfully completed.

Risk of injury if the device is operated outside the prescribed ambient conditions!

Use of the device outside the prescribed ambient conditions can result in failure to comply with tolerances, device failures, and injury to the patient.

⇒ Only operate the device within the prescribed ambient conditions (see chapter "Technical data").

⇒ Only operate the device indoors.

Notice

Dirt in the device can cause material damage!

Dirt entering the device during transport can damage the device.

⇒ Only transport the device with the system interface cover closed.

⇒ Transport the device in the corresponding transport bag.

2.2 General information

- Use of third-party products may lead to functional failures and restricted fitness of purpose. Biocompatibility may also be compromised. Please note that in these cases, any claim under warranty and liability will be void if neither the accessories nor original replacement parts recommended in the instructions for use are used.
- Repairs, servicing and maintenance should only be carried out by the manufacturer, or by a technician expressly authorized by the manufacturer.
- The device is subject to special precautions with regard to EMC (electromagnetic compatibility). Maintain the recommended safety distances (see ["10.2 Safety distances", page 45](#)) between the device and devices which emit high-frequency radiation (e.g. cellphones) to prevent malfunctions.
- Only connect up the approved SpO₂ sensors and the approved SpO₂ connection cable.
- Please note that no SpO₂ or pulse rate alarms occur in combination with WM 100 TD therapy devices.
- A functional test facility cannot be used to assess the accuracy of the SpO₂ values.
- When used with non-approved electrical medical devices or systems, increased electromagnetic emissions or reduced radio interference immunity may occur.
- The operator is responsible for ensuring the compatibility of the therapy device and all the connected components and accessories prior to the application with the patient. Only have modifications to the unit carried out by the manufacturer, or by a technician expressly authorized by the manufacturer.
- Please observe the chapter on hygienic preparation in order to avoid infection or bacterial contamination (see chapter "Hygienic preparation").
- Also observe the respective instructions for use for the therapy device, the components, and the accessories.

- When connecting the *VENTIremote alarm* or an internal hospital alarm system, ensure that any alarms and malfunctions can be perceived at all times in order to take the necessary measures.

2.3 Warnings in this document

Warnings are used to flag up safety-relevant information.

You will find a warning preceding any action that entails a hazard for persons or equipment.

Warnings consist of

- the warning symbol (pictogram),
- a signal word designating the hazard level,
- information about the hazard
- instructions for avoiding the hazard.

The warnings appear in three hazard levels, depending on the degree of danger:

Warning!

Designates an extremely dangerous situation. Failure to observe this warning may lead to serious, irreversible, or fatal injury.

Caution!

Designates a dangerous situation. Failure to observe this warning may lead to minor or moderately serious injury.

Notice!


Indicates a harmful situation. Failure to observe this warning may lead to damage to equipment.

Designates useful information relating to a particular action.

3 Product description


3.1 Module overview

3-1 prisma CHECK module


No.	Designation	Description
1	System interface cover	Covers the system interface if no additional module is connected.
2	SpO ₂ interface	Used to connect a pulse oximetry sensor.
3	System interface	Used for connection to the therapy device or other modules.
4	Locking hooks	Affix the module to the therapy device or other modules.
5	Module unlocking button	Allows connection of the module to the therapy device and removal of the module from the therapy device.
6	Remote alarm interface	Used to connect a remote alarm box or an internal hospital nurse call system. This function cannot be used with WM 100 TD therapy devices.
7	Mounting holes	Accept locking hooks of other modules.

3.2 Accessories

No.	Designation	Description
1	SpO ₂ connection cable	Connects the therapy device to the SpO ₂ sensor
2	SpO ₂ sensor	Used for connection to a patient's finger. Is available in different sizes.
3	Nurse call system connection cable	Connects the module to an internal hospital nurse call system
4	VENTiremote alarm connection cable	Connects the module to the remote alarm box VENTiremote alarm

3.3 Labels and symbols

3-2 Labels on the module


No.	Symbol	Description
Labels and symbols on the module		
1		Consult instructions for use
2		prisma CHECK module
4		Connection for pulse oximetry sensor
5		Connection for remote alarm (cannot be used with WM 100 TD therapy devices)
Labels and symbols on the type plate		
2	WM 100 MS	Type designation of prisma CHECK
	SN	Serial number
		Manufacturer and date of manufacture
		Type of protection against electric shock: Protection class II device

No.	Symbol	Description
2		Do not dispose of device in household waste
	IP21	IP protection class: Degree of protection against ingress of dangerous parts and solid foreign bodies. The device is protected against dripping water.
		Consult instructions for use
	CE 0197	CE mark (confirms that the product complies with the applicable European directives).

4 Preparation and operation

4.1 Connecting the module up to the therapy device

1. Open the system interface cover on the right side of the therapy device.


2. Keep the unlocking button on the module depressed.
3. Slide the module's locking hooks into the mounting holes on the therapy device until it clicks into place.

Result The module is connected.

4.2 Removing the module from the therapy device

1. Press the unlocking button on the module.


2. Remove the module sideways from the therapy device.
3. Close the system interface cover on the therapy device again.

Result The module is removed.

4.3 Connecting the pulse oximetry sensor

The pulse oximetry sensor and the SpO₂ connection cable are accessories and must be ordered separately.


1. Using the SpO₂ connection cable, connect the pulse oximetry sensor up to the SpO₂ interface.

CAUTION

Risk of injury if the sensor is fitted too tightly

An incorrectly fitted pulse oximetry sensor can occlude the blood supply.

- ⇒ Place the finger mark on the pulse oximetry sensor with the finger facing upwards.
- ⇒ Select the correct size of sensor.
- ⇒ Observe the instructions for use of the pulse oximetry sensor.

2. Place pulse oximetry sensor on one finger. Note: The sensor must be comfortably positioned and must not exert pressure on the finger.

Result The pulse oximetry sensor is connected. The display shows the current oxygen saturation and pulse rate values. 


The value under **SpO₂** displays the SpO₂ value in %. The value under **PR** displays the pulse rate in pulse beats per minute. If the values cannot be determined (e.g. because the sensor has not yet been connected or if lack of signal quality makes the values unreliable), - is displayed instead of the values.

4.4 Connecting the nurse call system and remote alarm box

Connect the prisma VENT to the nurse call system of a hospital or to the remote alarm box *VENTIremote alarm* via the remote alarm connection of the prisma CHECK module.

4.4.1 Connecting the remote alarm box *VENTIremote alarm*

Requirement The module is connected to a prisma VENT therapy device.

1. Connect the plug of the connection line to the remote alarm interface.
2. Connect the connection line to the remote alarm box *VENTIremote alarm* or the nurse call system. When doing so, observe the instructions for use for the remote alarm box or the nurse call system connection line.

4.4.2 Forwarding alarms

All high- and medium-priority alarms are switched in parallel to a remote alarm connection. The only low-priority alarm which is forwarded to the connection is the **No line voltage** alarm (if battery installed, otherwise high-priority). All other alarms are displayed only on the therapy device.

5 Hygienic preparation

5.1 General information

- Wear suitable protective equipment for disinfection work.
- Please refer to the instructions for use supplied with the disinfectant used.

5.2 Intervals

Interval	Action
Weekly	Clean the modules
As necessary, e.g., following infectious diseases or severe soiling	Disinfect the modules

5.3 Hygienic preparation of the module

NOTICE

Damage to the device caused by ingress of liquids!

Ingress of liquids may damage the therapy device, components, and accessories.

⇒ Do not immerse the therapy device, components, or accessories in liquids.

Requirement The module is disconnected from the therapy device (see "4.2 Removing the module from the therapy device", page 36).

1. Carry out hygienic preparation of the module as specified in the following table:

Part	Cleaning	Disinfection	Sterilization
Housing	Wipe with a damp cloth using water or mild soap	Wipe disinfection (Recommendation: terralin® protect or perform advanced Alcohol EP)	Not permitted
Connection cable			

2. Allow the module housing and in particular the electrical interfaces to dry well.

Result The module is hygienically prepared.

6 Function check

6.1 Intervals

Carry out a function check at regular intervals:

- After each hygienic preparation
- After each repair
- At least every 6 months

6.2 Performing a function check

1. Inspect the housing for cracks, damage, and heavy soiling.
2. Connect the module up to the therapy device (see "4.1 Connecting the module up to the therapy device", page 35).
3. Check whether the  symbol is shown on the therapy device display.
4. Connect the pulse oximetry sensor (see "4.3 Connecting the pulse oximetry sensor", page 37).
5. Check whether the current oxygen saturation and pulse rate values are shown on the display.
6. If the module does not work properly or shows signs of damage: Contact your authorized dealer.

Result The function check is complete.

7 Faults

If you are unable to remedy faults with the aid of the table, or in the event of unexpected operation or an incident, contact your specialist dealer to have the device repaired. To avoid serious damage, do not continue using the device.

Fault	Cause	Remedy
 does not appear on the therapy device's display.	Contact problem or defective module.	Remove and then reconnect the module (see 4.2, p.36), (see 4.1, p.35). If the fault persists: Have the module replaced.
No values for SpO ₂ and PR are displayed.	Contact problem or defective SpO ₂ sensor.	Disconnect the sensor from the SpO ₂ connection cable or the SpO ₂ connection cable from the module and reconnect. If the fault persists: Have the sensor replaced.
	The SpO ₂ sensor is poorly positioned	Check the position and fitting of the SpO ₂ sensor. Select another sensor size if necessary.
	Nail polish or contaminated finger	Check finger. If necessary: Clean finger.

8 Maintenance

The module is designed to have a useful service life of 6 years.

If the module is used as intended in accordance with these instructions for use, it does not require any maintenance.

If you identify faulty parts during the function check ([see "6 Function check", page 41](#)), contact your authorized dealer.


9 Storage and disposal

9.1 Storage

Store the device under the specified ambient conditions (see chapter "Technical data").

9.2 Disposal

9.2.1 Electronic waste

Do not dispose of the product in the household waste. Consult an authorized, certified electronic waste recycling company for proper disposal. You can find out their address from your environmental officer or from your local council.

The device packaging (cardboard box and inserts) can be disposed of as waste paper.

10 Appendix

10.1 Technical data

Specification	prisma CHECK
Product class according to 93/42/EEC	IIa
Dimensions W x H x D in cm	4.3 x 12.8 x 8
Weight in kg	0.185
Temperature range Operation Storage	+5°C to +40°C -25°C to +70°C
Permissible humidity during operation and storage	Rel. humidity 10% to 90%, non-condensing
Air pressure range	700 hPa to 1060 hPa, corresponds to a height of approx. 3000 m above sea level
Power supply	12 V DC
Electrical output	1.2 VA*
Insulation of remote alarm interface	1.0 kV
Switching capacity of remote alarm interface	60V DC/2A; 42V AC/2A
SpO ₂ display Pulse display Update period for SpO ₂ and pulse	0% to 100% 0 bpm to 300 bpm 0.5 s
Applied part	SpO ₂ sensor
Classification as per DIN EN 60601-11: Type of protection against elec. shock Degree of protection against elec. shock Protection against harmful ingress of water and solid particles	Protection class II Type BF IP21
Classification as per DIN EN 60601-1: Operating mode	Continuous operation
Electromagnetic compatibility (EMC) as per EN 60601-1-2 Radio interference suppression Radio interference immunity	Electrical medical devices may only be installed and commissioned in a defined electromagnetic environment with regard to emission and radio interference immunity. More information, including test parameters and limit values, can be obtained from the manufacturer if required. EN55011 B IEC 61000-4 Parts 2 to 6, Part 11, Part 8 IEC 61000-3 Parts 2 and 3

*Data reached only in combination with the approved therapy device.

The right to make design modifications is reserved.

All flow values are determined under STPD conditions.

10.2 Safety distances

Recommended safety distances between portable and mobile HF telecommunication devices (e.g. cell phones) and the device in order to prevent malfunctions.			
Nominal capacity of HF device in W	Safety distance depending on transmission frequency in m		
	150 kHz - 80 MHz	80 MHz - 800 MHz	800 MHz - 2,5 GHz
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.69	3.69	7.38
100	11.67	11.67	23.33

10.3 Scope of supply

prisma CHECK module	WM 29390
Part	Article number
prisma CHECK basic device	WM 29391
Instructions for use	WM 68190

10.4 Accessories

Accessories can be ordered separately, if required. A current list of scopes of supply can be ordered on the website of the manufacturer or through your specialist dealer.

10.5 Spare parts

Replacement parts can be ordered separately, if required. A current list of accessories and replacement parts can be ordered on the internet site of the manufacturer or through your authorized specialist dealer.

10.6 Warranty

Starting from the date of purchase, Löwenstein Medical offers the customer a limited manufacturer's warranty on a new original Löwenstein Medical product or replacement parts installed by Löwenstein Medical in accordance with applicable warranty terms and conditions for the particular product and the warranty periods listed below. The warranty conditions are available on the website of the manufacturer. We can also send you the warranty conditions on request.

In the event of a claim under warranty, please contact your specialist dealer.

Product	Warranty periods
Löwenstein Medical devices, incl. accessories (excluding: masks), for sleep diagnostics, home ventilation, oxygen therapy, and emergency medicine	2 years
Masks, incl. accessories, batteries (unless otherwise stated in the technical documentation), sensors, hose systems	6 months
Disposable products	None

10.7 Declaration of conformity

Löwenstein Medical Technology GmbH + Co. KG, Kronsaaßweg 40, 22525 Hamburg, Germany, the manufacturer of the devices described in these Instructions for Use, hereby declares that the product complies with the respective regulations of Medical Devices Directive 93/42/EEC. The unabridged text of the Declaration of Conformity can be found on the manufacturer's website.

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